

# Your Natural Choice for the Cardiac Cath Lab

Enhance your ability to provide angiographic care with an injector that delivers versatility, control, interoperability<sup>1</sup> – and the trusted support of a pioneering leader in the field of angiography.

The Bayer MEDRAD® Avanta injector is an established contrast injection system to lead your clinical success in angiographic diagnostic and interventional procedures.

Surveyed interventional cardiologists are more likely to recommend the Bayer MEDRAD® Avanta injector than the leading competitor injector.<sup>2</sup>



## Your Versatile IC Injector

The Bayer MEDRAD® Avanta offers all-around flexibility.

Five patients  
with one set-up



The multi-patient tube, contrast, syringe, and hand controller can be used for up to five consecutive patients – saving time between patients.

Fluid management  
system



For precise delivery of contrast and saline throughout your diagnostic and interventional procedures.


High customer  
satisfaction



Surveyed interventional cardiologists are highly satisfied with Avanta for performing PCI (percutaneous coronary interventions) and TAVI/EVAR procedures.<sup>2</sup>

Clear Direction.  From Diagnosis to Care.

**MEDRAD® Avanta**  
Fluid Management Injection System



## Precision and Control

Precision control of fluid delivery, together with inline air detection, gives you confidence for your interventional procedures.

- The dual-line tube enables real-time control of contrast media and saline – before or after injection. The hand controller helps achieve an accurate flow rate, as well as giving you the option of manually controlled injections. In addition, the MEDRAD® Avanta delivers a precise contrast media volume, as well as bolus sharpness for controlled contrast delivery.<sup>1</sup>



## Easier Workflow Integration

Because each suite is unique, with MEDRAD® Avanta, you can choose the configuration that optimises your workflow.

- The pedestal option offers mobility of the injector, with a large handle to help you easily move the injector.
- With the table mount option, the injector moves with the table, offering an integrated workflow.
- Sterile disposables provide the precise fluid control necessary for coronary imaging, as well as a single set-up that allows for up to five patients, further support the needs of your cardiac cath lab.

## From a Trusted Pioneering Leader in The Field of Angiography

Bayer has over 50 years of experience in the field of power injection and contrast media.

- On-site clinical training to help you perform procedures with confidence from day one
- On-site device service and remote technical support
- Our quality care comes in customised solutions

1. Bayer. Avanta Operation Manual. 2021.
2. Bayer, internal data on file. Blinded Customer Perception Survey among 160 IR + IC, in US, DE, IT, and JP, 2021.file. Blinded Customer Perception Survey among 160 IR + IC, in US, DE, IT, and JP, 2021.

Clear Direction.  From Diagnosis to Care.

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorized Bayer representative for the most current information.

The patient data that appears in this document is actual health information but all personal identifiers have been removed or otherwise anonymized. No personally identifiable information is shown.

Bayer, the Bayer Cross, MEDRAD Avanta, MEDRAD, Avanta, Ultravist, VirtualCare, DirectCare, PartnerCare, and SelectCare are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

© 2024 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.



More information on  
[radiology.bayer.com.au](http://radiology.bayer.com.au)



**Manufacturer**

Bayer Medical Care Inc.  
1 Bayer Drive Indianola, PA  
15051-0780 U.S.A.  
Phone: +1-412-767-2400  
+1-800-633-7231  
Fax: +1-412-676-4120

**Australian Sponsor**

Imaxeon Pty Ltd  
Unit 1, 38 – 46 South Street  
Rydalmere, NSW 2116 Australia  
Phone: 02 8845 4999  
Fax: 02 8845 4936  
Customer Service: 1800 633 723