



Workflow Solutions //Injector Management

Because Best Practice Is Best Shared

Clear Direction.  From Diagnosis to Care.

Share Protocols at the Press of a Button

Injection protocols have a direct influence on the result of contrast-enhanced exams. With a larger number of injection systems, it may become increasingly difficult to keep all devices harmonised. This may impact the quality and consistency of contrast procedures.

//Injector Management allows for central administration of protocols, catheter, and contrast libraries. With the radiology team in control of their best practices, the department can set its own standards and operate efficiently and successfully.






Set Your Standards

- Sharing best practices is now just a click away to help ensure adherence to standards in the clinical routine.



Full Transparency

See all libraries on connected devices in a single view for convenient access and a concise overview.



Group Multiple Devices

Multiple devices can be grouped together to share identical libraries and to simplify management.



Define New Profiles

Create and manage profiles by either duplicating your best injector's profile or defining custom libraries.



Learn One, Know All

- Grouped injectors share common protocols and names.



Enhanced Confidence

Confidently use familiar protocols within a group.



Easy Learning

Identical injector protocol nomenclature across scanner suites means easy learning, allowing for simplified staff rotation.



Full Flexibility

All injectors retain the possibility to modify settings at any time.



Better Together

➤ Networking injectors benefits the whole team.



Suggest New Protocols

Best practice can still evolve so //Injector Management makes it easy to suggest new protocols for inclusion in the library.



Centrally Stored Backup

Exchanging injectors for service is faster and simpler with a centrally stored backup.



Patient Priors

When combined with the available //Smart Protocols software, patient priors from connected devices are selectable, e.g. for follow-ups.

//Injector Management

Share injector protocols through a centralised server to make them available to connected injectors.


Central Server Unit

- Central administration of protocol, fluid, catheter, and injector sites libraries
- Review and approve protocols
- Centralised backup
- Group up to 25 injectors together for consistent settings



Suggest Protocols

- Users can suggest protocols for inclusion in the library



Injector Groups

Identical libraries in injector groups for

- Protocols
- Fluid
- Catheter
- Injection Sites



Clear Direction. From Diagnosis to Care.

Bayer reserves the right to modify the specifications and features described herein or to discontinue any product or service identified in this publication at any time without prior notice or obligation. Please contact your authorised Bayer representative for the most current information.

All patient data that appears in this document are fictitious. No actual patient information is shown.

Bayer, the Bayer Cross, MEDRAD Stellant, MEDRAD MRXperion, MEDRAD Stellant FLEX, MEDRAD, Stellant, MRXperion, Stellant FLEX are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

© 2022 Bayer. This material may not be reproduced, displayed, modified or distributed without the express prior written consent of Bayer.



Bayer Healthcare LLC
100 Bayer Boulevard
P.O. Box 915
Whippany, NJ 07981
U.S.A

More information on
radiology.bayer.com.au



Manufacturer Stellant FLEX
& MRXperion
Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.

Phone: +1-412-767-2400
+1-800-633-7231
Fax: +1-412-676-4120

Manufacturer Centargo
Imaxeon Pty Ltd
Unit 1, 38-46 South St
Rydalmere, NSW 2116
Australia

Phone: 02 8845 4999
Fax: 02 8845 4998
Customer Service: 1800 633 723

Australian Sponsor
Imaxeon Pty Ltd
Unit 1, 38-46 South St
Rydalmere, NSW 2116
Australia

Phone: 02 8845 4999
Customer Service: 1800 633 723